

ASTM LEVEL 1 ISOLATION MASK


FDA Certified Isolation-Grade Medical 3-Ply Mask



3-Layer protection designed to effectively block out dust, bacteria, droplets, and air pollution.

ASTM International sets the preferred international standard in healthcare for best practices inclusive of testing and requirements for performance of materials.

The FDA, who oversees the sale and marketing of all surgical masks, recommends that manufacturers demonstrate mask performance in 4 areas: Fluid Resistance, Differential Pressure, Filter Efficiency, Flammability



Built-in adjustable metal nose-bridge strip, to better fit and seal your face

<p>ASTM LEVEL 1</p> <p>Low Fluid Resistance 80 mmHg Filtration Efficiency BFE ≥ 95% PFE ≥ 95% @ 0.1 micron < 4.0 mm H₂O/cm² Breathability - Delta P Flame Spread Class 1</p>		<p>Ideal for procedures where low amounts of fluid, spray and/or aerosols are produced.</p> <p>Meets EN14683 Rating – Type II Standard.</p> <p>Pictured: Isofluid® Earloop with Secure Fit™ Mask Technology</p>	
<p>LOW PERFORMANCE</p> <p>Surgical Molded Utility Mask Physical Barrier Only No LEVEL Performance Level ** Filtration Efficiency N/A</p> <p>**Unless mask manufacturer certifies mask meets ASTM performance Level 1</p>		<p>Ideal as a comfortable substitute for earloop face masks, this mask is a simple physical barrier for exams and visitations or for dry, short procedures that do not produce fluid, spray or aerosols.</p> <p>Pictured: Surgical Molded</p>	
<p>MINIMUM PERFORMANCE</p> <p>Utility Mask (Tissue/Tissue) Physical Barrier Only No LEVEL Performance Level Filtration Efficiency N/A</p>		<p>Ideal as a simple physical barrier for exams and visitations or for dry, short procedures that do not produce fluid, spray or aerosols.</p> <p>Pictured: Isolite® Earloop</p>	

Follow CDC Guidelines: Do not treat active TB patients except in approved facilities, meeting all health department, CDC and OSHA standards, in the context of a complete respiratory protection program. Wash hands after touching any contaminated surfaces. Do not touch outside of the mask with wet or contaminated gloves or hands. Such contamination may compromise mask barrier asepsis by encouraging migration or “wicking” of microbes through the mask. Mickey’s Medclean, can make no warranties or representations, either expressed or implied, that these products will fully protect the user from exposure to blood or bodily fluids or risk of contracting a virus or infectious diseases.



Delivery dates are best estimates.

Pricing is valid on day of quote, and may fluctuate due to very limited supplies and extreme surges in air freight.